

Development of a Multi-Peptide Immunocapture – LC-MS/MS Assay for the Quantitation of a **PEGylated Therapeutic in Rat and Dog Plasma**

Kevork Mekhssian¹, Jean-François Dupuis¹, Olivier Didur², William Ruddock² and Anahita Keyhani¹ ¹Altasciences, Laval, Quebec, Canada; ²ITR Laboratories, Montreal, Quebec, Canada

OVERVIEW

PURPOSE

To develop a multi-peptide assay based on LC-MS/MS for the quantitative analysis of a PEGylated therapeutic protein to support preclinical studies.

METHOD

Pegvisomant was extracted from plasma using direct digestion or immunocapture followed by on-bead tryptic digestion. Two surrogate peptides, N-term LFDNAMLR and C-term VSTFLR peptides were used for quantitation.

RESULTS

Both direct digestion and immunocapture with anti-PEG showed good linearity, precision and accuracy. However, higher sensitivity and specificity were obtained with the immunoaffinity approach. Good agreement of Pegvisomant concentration values between VSTFLR and LFDNAMLR-based quantitation were observed in rat and dog PK samples.

INTRODUCTION

While advances in HRMS have facilitated the analysis of the intact form of proteins, the preferred strategy for the quantitative bioanalysis of large molecules remains the surrogate peptide approach. However, metabolic protein cleavage followed by potentially divergent degradation pathways could introduce a bias in the ultimate determination of intact analyte concentration. To address this challenge, an immunocapture liquid chromatography-mass spectrometry (IC-LC/MS) analysis strategy was developed and optimized for quantitative analysis of a PEGylated therapeutic. The strategy consisted in selecting two surrogate peptides from opposite terms. This dual-peptide quantitation approach was put to test by generating single dose pharmacokinetic profiles of Pegvisomant (1) in both rat and dog models.






Figure 1. Amino acid sequence of Pegvisomant protein

METHODS

PEGVISOMANT EXTRACTION PROCEDURE



IMMUNOCAPTURE



LC-MS/MS DETECTION

Table 1. Pegvisomant method MRM transitions.

LFDNAM VSTFL VSTFLR[^]

PHARMACOKINETIC STUDY

 Four Sprague Dawley rats and one Beagle dog Subcutaneous dose at 2 mg/kg • Collection at 0,1, 8, 24, 48, 72, 96 and 120 hours

Biotinylated α PEG and α hGH coupled to Strep magnetic beads Magnetic sample processing using KingFisher FlexTM

Agilent Zorbax SB300 C₁₈ column (50 x 2.1mm, 3.5 μ m) • Gradient elution with 0.1% CH₃COOH in H₂O and ACN SCIEX Triple Quad API 5000 operated in (+) ESI-MRM

	Q1 m/z	Q3 m/z	CE (eV)	Dwell Time (ms)
R	490.2 (2+)	719.3 (y6)	20	50
R	361.7 (2+)	536.3 (y4)	20	50
(IS)	366.7 (2+)	546.5 (y4)	20	50

RESULTS

DIRECT DIGESTION vs. IMMUNOCAPTURE LC-MS

Both direct digestion and immunocapture LC-MS approaches were evaluated. Although the direct digestion method is fast and low-cost, it is less sensitive and prone to interferences from the matrix.

Direct digestion of Pegvisomant in plasma after reduction-alkylation shows maximum MS response for both peptides after 20 min incubation with porcine trypsin at 37 °C (Figure 3).






Figure 3. Direct digestion kinetics of Pegvisomant using porcine trypsin

For the immunocapture approach, higher extraction efficiency is obtained using anti-PEG antibodies compared to anti-hGH (Figure 4).



various anti-PEG and anti-hGH antibodies

LFDNAMLR



Significant increase in sensitivity is achieved using the immunocapture PK APPLICATION approach with anti-PEG mAb for Pegvisomant extraction of only 10 μ L Concentrations of Pegvisomant in PK samples (Table 2) generated using of rat plasma compared to 25 μ L using the direct digestion approach either N-term or C-term peptide show good agreement, confirming the (Figure 5). structural integrity of the PEGylated protein (Figure 7).




Figure 5. Chromatograms of extracted LLOQ (25 ng/mL) for Pegvisomant in Rat plasma using direct digestion (top) and ar PEG immunocapture (bottom)

Samples extracted using immunocapture show excellent linea precision in the analytical range (25 - 10,000 ng/mL) (Figure 6).




Figure 6. Calibration curves of Pegvisomant from 25 ng/mL to 10,000 ng/mL in rat plasma using the immunocapture approach

LEARN MORE ABOUT OUR BIOANALYSIS AND RESEARCH SERVICES ()

Altasciences, 575 Armand-Frappier, Laval (Montréal), Québec, Canada contact@altasciences.com altasciences.com




Figure 7. PK profiles of Pegvisomant in Sprague Dawley rat (n=4) (left) and Beagle dog (n=1) (right) following single SC dose at 2.0 mg/kg

	Conc. (ng/mL)	Time (hr)						
1.8		1	8	24	48	72	96	120
	VSTFLR _Rat	388	1663	5914	6764	4768	3108	2009
anti-	LFDNAMLR_Rat	418	1592	5379	5999	4608	2949	2093
	Difference (%)	+7.4	-4.4	-9.5	-12.0	-3.4	-5.2	+4.1
earity and	VSTFLR_Dog	999	5576	6843	7811	8398	8601	9923
).	LFDNAMLR_Dog	1019	5726	7067	8466	8665	9238	10177
	Difference (%)	+2.0	+2.7	+3.2	+8.0	+3.1	+7.1	+2.5

Table 2. Pegvisomant Conc. following SC 2 mg/kg dose in rat and dog

CONCLUSION

The strategic selection of multiple surrogate peptides for protein quantitation can help address some of the issues facing the bottom-up bioanalysis of proteins by LC-MS and in this case, ensure the structural integrity of Pegvisomant in different biological matrices. Furthermore, this study highlights the advantages of immunocapture LC-MS compared to the direct digestion approach for the quantitation of biotherapeutics.

REFERENCES

Trainer et al. Treatment of acromegaly with the growth hormone-receptor

agonist pegvisomant. N Engl J Med 2000; 342:1171-7.