

# Development of a Hybridization LC-MS Methodology for Quantitation of siRNA in Cynomolgus Monkey


Jean-François Dupuis<sup>1</sup>, Long Yuan<sup>2</sup>, Jean-Nicholas Mess<sup>1</sup> and Kevork Mekhssian<sup>1</sup>

<sup>1</sup>Altasciences, Laval, QC, Canada; <sup>2</sup>Biogen, Cambridge, MA, USA

[Click here to listen to the recorded poster presentation](#)

## INTRODUCTION


Small interfering RNA (siRNAs) are short, double-stranded oligonucleotides that can regulate gene expression by inducing mRNA degradation (**Figure 1**). Quantitative bioanalysis in various biological matrices is critical to understanding the PK/PD parameters of therapeutic oligonucleotide drug candidates. Hybridization LC-MS/MS has been recently established as a successful methodology for quantifying single-stranded oligonucleotides in plasma. This approach relies on capturing the oligonucleotide of interest with a full-length complementary DNA probe, followed by analysis with LC-MS/MS (**Figure 2**). Applying the same workflow to double-stranded siRNA presents additional challenges, as the probe competes with the sense-strand of the siRNA. This study presents the first adaptation of the hybridization LC-MS/MS workflow for quantifying a double-stranded oligonucleotide (siRNA1).




**Figure 1.** siRNA mechanism of action (Image adapted from Petrova N.S *et al.*, 2013)

## METHODS

### HYBRIDIZATION LC-MS/MS



**Figure 2.** Hybridization LC-MS/MS workflow for the quantitation of siRNA in biological matrices




**Figure 3.** Deoxyribonucleic acid (DNA) and peptide nucleic acid (PNA) structures

**Table 1.** Analyte and Internal Standard Parameters

Oligo	MW (kDa)	Length (bp)	Chemistry	MRM
AS1 (siRNA1 antisense strand)	6.9	21	Mixed backbone with 2'-OMe and 2'-F modification	629 / 95
Internal Standard	7.9	20	Uniform MOE with PS backbone	879 / 95


## RESULTS

In the hybridization step, a full-length complementary DNA-based probe is typically employed to capture single-stranded oligonucleotides because of its high affinity and specificity. This DNA-based probe achieved near-complete recovery of the antisense strand (AS1) of siRNA1, independent of temperature and salt concentration (**Figure 4, left**). However, under the same conditions, recovery of the double-stranded siRNA was low at lower temperatures and high salt concentrations due to the competition from the siRNA1 sense-strand (**Figure 4, right**).



**Figure 4.** Capture efficiency (% Recovery) of AS1 (left) and siRNA1 (right) using a DNA-based probe at different temperature and salt concentrations

The recovery of siRNA1 utilizing PNA and DNA-based probes was compared. The PNA capture probe offers a key advantage due to its lack of a negatively charged deoxyribose phosphate backbone, which eliminates electrostatic repulsion with the siRNA1 antisense strand. Moreover, the higher melting point of the PNA:RNA complex allows for more efficient annealing at elevated temperatures. Under conditions of 65 °C and 10 mM NaCl (**Figure 5**), the PNA probe achieved a 90% recovery compared to 55% with the DNA probe. The PNA-probe conjugated beads were then incorporated into the extraction workflow, and the assay was evaluated for precision and accuracy (**Table 2**).



**Figure 5.** Capture efficiency (% Recovery) of siRNA1 using PNA- versus DNA-based probe at increasing temperature

**Table 2.** Precision and Accuracy of siRNA1 in Cyno Plasma Extracted Using the PNA Probe

QC Sample	LLOQ QC 2.00 ng/mL	LQC 6.00 ng/mL	MQC 500 ng/mL	HQC 750 ng/mL
1	1.98	5.87	527.11	681.34
2	1.94	6.50	441.75	674.75
3	2.05	5.91	464.04	756.63
4	2.29	5.56	469.18	753.45
Mean	2.07	5.96	475.52	716.54
S.D.	0.16	0.39	36.40	44.55
n	4	4	4	4
% C.V.	7.6	6.6	7.7	6.2
% Nominal	103.3	99.3	95.1	95.5

## CONCLUSION

An analytical method for the quantitation of double-stranded siRNA1 using a PNA-based probe was successfully developed. The method is accurate and precise across a range of 2.00 to 1000.00 ng/mL in cyno plasma. With minor adjustments, the assay was subsequently qualified in CSF and tissues and used to support a nonclinical study in cynomolgus monkeys.

## REFERENCE

Yuan, L. Dupuis, J-F, Mekhssian, K. A Novel Hybridization LC-MS/MS Methodology for Quantification of siRNA in Plasma, CSF and Tissues Samples *Molecules* 2023